**BRESSER®**

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

**GB Instruction Manual**

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
**BRESSER Science MPO-401**  
**Transmission-type Polarization Microscope**

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### **CAUTION!**

To work with this microscope, sharp and pointed aids are being used. Please take care that this microscope and its accessories are stored at a place out of reach of children. Let children only work with this microscope under an adult's supervision! Keep packing material (plastic bags etc.) away from children!

### **Note!**

If you have any complaints or queries please first contact your national service centre by telephone. The address is included in these instructions.

### **Only for EU countries**

**Do not dispose of electric equipment together with household waste material!**

In observance of European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE) and its implementation in accordance with national law, electric equipment that have reached the end of their life must be collected separately and returned to an environmentally compatible recycling facility.

Discharged batteries and damaged re-chargable batteries must be disposed of at special battery collection points. Information is available from your local disposal agent or local authority regarding the disposal of devices or batteries manufactured after the 01.06.2006

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## General Information

### About this Instruction Manual

Please read the safety instructions in this manual carefully. To avoid damage to the unit and the risk of injury, only use this product as described in the manual.

Keep the instruction manual handy so that you can easily look up information on all the functions.

#### **DANGER!**

You will find this symbol in front of every section of text which deals with the risk of severe injury or even death in the event of improper use.

#### **CAREFUL!**

You will find this symbol before every section of text that deals with the risk of minor to severe injuries resulting from improper use.

#### **NOTE!**

You will find this symbol in front of every section of text which deals with the risk of damage to property or the environment.

### Intended Use

This product is intended only for private use. It was developed for the magnified display of things in nature.

### General Warning

#### **DANGER!**

Tools with sharp edges and points are often used when working with this device. For this reason, store this device and all accessories and tools in a location that is out of the reach of children. There is a RISK OF INJURY!

#### **DANGER!**

This device contains electronic components which operate via a power source (power supply and/or batteries). Do not leave children unattended when using the device. Only use the device as described in the manual, otherwise you run the RISK of an ELECTRIC SHOCK.

#### **DANGER!**

Do not expose the device to high temperatures. Use only the supplied power adapter. Do not short circuit the device or batteries, or throw them into a fire. Excessive heat or improper handling could trigger a short circuit, fires or even explosions.

#### **DANGER!**

Never bend, pinch, pull or run over the power and connecting cables or extensions/adapters. Protect the cables from sharp edges and heat. Before operation, check the device, cables and connections for damage. Never use a damaged unit or a unit with damaged power cables, etc. Damaged parts must be exchanged by an authorized service centre immediately.

Position your device so that it can be disconnected from the power supply at any time. The wall socket should always be located near the device and be easily accessible, since the plug on the power cord serves as a disconnecting device for the power supply.

#### **DANGER!**

Children should only use the device under adult supervision. Keep packaging material (plastic bags, rubber bands, etc.) out of the reach of children! There is a RISK OF CHOKING!

#### **CAREFUL!**

Children must not have access to the included chemicals and liquids. Do not drink the chemicals. Wash hands thoroughly with running water after use. In the event of contact with the eyes or mouth, rinse thoroughly with water. In the event of pain, contact a doctor immediately and take the substances with you.

#### **NOTE!**

Do not disassemble the device. In the event of a defect, please contact the Service Centre in your country (see the warranty card).

Do not expose the device to temperatures above 40°C.

### NOTES on cleaning

Separate the device from the power supply before cleaning (remove plug).

Only use a dry cloth to clean the exterior of the device.

#### **NOTE!**

Do not use any cleaning fluid to avoid damaging the electronics.

Clean the lenses (eyepiece and/or lens) with a soft and lint-free cloth only (e.g. microfibre).

#### **NOTE!**

Do not apply excess pressure to the cloth so as to avoid scratching the lenses.

To remove more stubborn dirt, moisten the cleaning cloth with an eyeglass-cleaning solution and wipe the lenses gently.

Protect the device from dust and moisture! After use - in particular in situations of high humidity - let the device acclimatize for a short period of time, so that the residual moisture can dissipate. Remove the dust cover and store it in the included bag.

### Disposal

Dispose of the packaging materials properly, according to their type (paper, cardboard, etc). Contact your local waste disposal service or environmental authority for information on the proper disposal.

Do not dispose of electronic devices in the household garbage! As per the Directive 2002/96/EC of the European Parliament on waste electrical and electronic equipment and its adaptation into German law, used electronic devices must be collected separately and recycled in an environmentally friendly manner. Empty old batteries must be disposed of at battery collection points by the consumer. You can find out more information about the disposal of devices or batteries produced after 01.06.2006 from your local waste disposal service or environmental authority.

Please take the current legal regulations into account when disposing of your device. You can get more information on the proper disposal from your local waste disposal service or environmental authority.

Please disconnect the plug from main socket before changing the fuse, removing the bottom on the rear plate or opening the lamp door to change the lamp.

The power cord provided with the equipment has a grounded plug. Always use the power cord with a properly grounded mains socket.

Position your device so that it can be disconnected from the power supply at any time. The wall socket should always be located near the device and

be easily accessible, since the plug on the power cord serves as a disconnecting device for the power supply.

Do not expose the instrument to high temperatures or humidity. Avoid using the instrument in extremely dusty locations.

Operating temperature: 5 °C to 35 °C.

Operating humidity: 20% to 80% (at 25 °C).

## Table of Contents

- I. Components
- II. Technical Data
- III. Setup Instructions
- IV. Observing Operating
- V. Maintenance, EEC Conformity Explanation, Warranty

### I. Components

Fig. 1

- 1 Camera port
- 2 Microscope head set screw
- 3 Bertrand lens centering screw
- 4 Holding clamp for slide
- 5 Condenser (with aperture diaphragm)
- 6 Polarizer
- 7 Field diaphragm (with filter seat)
- 8 Brightness control knob
- 9 Power switch
- 10 Fine focusing knob
- 11 Knob for tensional adjustment of focusing
- 12 Condenser up-down knob
- 13 Stage centering screw

Fig. 2

- 1 Eyepiece
- 2 Trinocular head
- 3 Bertrand lens
- 4 Analyzer
- 5 Intermediate attachment set screw
- 6 Limit knob (up stop)
- 7 Coarse focusing knob
- 8 Power cord
- 9 Mains in
- 10 Objective
- 11 Nosepiece
- 12 Compensator
- 13 Intermediate attachment

### III. SETUP INSTRUCTIONS

1. Remove all parts from their packing materials and retain the packaging in the event you need to transport the product.
2. Loosen the optical head set screw and turn the trinocular head to the face of the instrument. Secure the optical head with the screw.
3. You may install the Bertrand lens, the analyzer and compensator to the intermediate attachment.

## II. TECHNICAL DATA

### Eyepieces

Eyepiece type	Magnification	Field of view $\varnothing$ / mm	Focal length $f$ / mm	Remarks
Wide field eyepiece (WF)	10x	18	25	
Wide field eyepiece (WF) with crosshair	10x	18		optional

### Objectives

Objective type	Magnification	Numerical aperture (NA)	Remarks
Polarization-objectives (POL)	4x	0.1	
	10x	0.25	
	40x	0.65	
	100x	1.25 (oil)	
Planachromatic objective (without covering glass usable)	4x	0.1	optional

### Total magnification

Objective:	4x	10x	40x	100x
Total Magnification				
Eyepiece:				
10x	40x	100x	400x	1000x

- Power supply: AC 230 V / 50 Hz.
- Illumination: Halogen lamp (6 V, 20 W) with dimmer.
- Polarizer and analyzer.


- Rotable stage.
- Fine focusing knob: Minimum division 0.002 mm.
- Interpupillary distance adjustment range: 53 mm - 75 mm.

Reservation of technical alterations

4. Install the 40x objective to the nosepiece, into the opening without centering screw. Then install the other objectives one by one.
5. Remove the dust caps from the eyepiece holders on the optical head and insert both eyepieces into their holders.
6. Connect the power cord to a suitable power supply.

#### IV. OBSERVING OPERATIONS

##### 1. Illumination Controls (Fig. 3)

- 1 Power switch
- 2 Brightness control knob

1. The power switch for the illuminator and the brightness control are located on the base. The electrical system is fuse-protected and the fuse holder is located on the mains in.
2. Turn on the light with the power switch. If the light does not appear to be ON, check the brightness control to see if it's on a sufficiently low setting.

##### 2. Focusing Controls (Fig. 4)

- 1 Coarse focusing knob
- 2 Fine focusing knob
- 3 Knob for tensional adjustment of focusing
- 4 Limit knob (up stop)

1. **Focusing adjustment:** Focusing adjustment is accomplished by using the large coarse adjustment knobs (1) located comfortably on each side of the frame. Fine adjustment is accomplished using the smaller knobs (2) located on the same focusing shaft. This coaxial arrangement allows for easy, precise adjustment without drift or discomfort.
2. **Focus control:** Turning either of the focusing knobs will raise or lower the stage. The smallest graduation on the index scale of the fine focusing knob (2) is 2  $\mu\text{m}$  of vertical.
3. **Focusing tension adjustment:** The tension of the coarse focusing knob is adjustable and preset at the factory for ease of use. If you wish to adjust the coarse focusing tension, first locate the ringed tension adjustment knob (3).

It is located between the frame and the left coarse focusing knob. Turning the ring toward the front of the microscope increases the tension, and toward the rear of the microscope loosens it. The tension is too high if you experience physical discomfort.

4. **Pre-focusing or focus stop control:** The use of this feature will insure that the shorter working distance objectives don't contact the stage or specimen when using the microscope. Its use also simplifies focusing. After focusing on the specimen with the coarse adjustment by a low magnification eyepiece, turning up and fastening of the lever (4) will set an upper limit on the coarse focusing movement. After changing specimens or objectives, focusing is easily accomplished by rotating the coarse focusing knob to reach the pre-focused position. Then make fine adjustments with the fine focusing knob. Focusing movement with the fine focusing knob isn't affected by using the pre-focusing lever.

##### 3. Interpupillary and Diopter Adjustments (Fig. 5)


- 1 Binocular tube
- 2 Diopter adjustment ring

1. **Interpupillary distance adjustment:** Proper interpupillary distance, or the distance between the eyepieces, is crucial to the comfort of the user. Adjusting the interpupillary distance is accomplished through a „folding“ action of the binocular tube (1), allowing for quick and easy adjustment.
2. **Diopter adjustment:** Proper correction for the individual vision is accomplished via the diopter adjustment located at the left eyepiece. Using the objective with 40x magnification, close your left eye and bring an image into focus with your right eye only. Once the image is well focused, observe with the left eye and make fine adjustments with the diopter adjustment ring (2) to correct for your vision.


#### 4. Rotatable Stage Controls (Fig. 6, A, B, C)

- 1 Stage centering screw
- 2 40x objective
- 3 Objective centering screw
- 4 Stage set screw

1. Put a specimen on the stage, observe it using the crosshair eyepiece and the 40x objective (2).


- 1 Target
- 2 Center of view
- 3 Center of circle
- 4 Circle

2. Search one target in the field of view, adjust it to the center of view (see Fig. A).
3. Rotate the stage. If the stage is out of center, the target will encircle the center of an eccentric circle (see Fig. B).
4. In order to center the center of the circle in the the field of view, use the two stage centering screws (1). The stage centering adjustment finished, see figure C.
5. If other objectives are out of center, you may adjust them using the objective centering screws (3) driven by the Allen-wrenches.
6. Tightening the stage set screw (4) can fix the rotatable stage.

#### 5. Condenser Alignment Controls (Fig. 7, A, B, C)


1. Condenser Components: The components of the Abbe condenser (Fig. 1, No. 5) are located below the stage. There are its mounting ring with holding screw (4) and two large centering screws (5) at either side. At the left, there is an up-down knob (1) for raising and lowering the condenser.

- 1 Condenser up-down knob
- 2 Auxiliary lens of condenser
- 3 Aperture adjustment lever
- 4 Condenser holding screw
- 5 Condenser centering screw
- 6 Polarizer
- 7 Field diaphragm
- 8 Filter seat

2) Condenser Alignment: This procedure is done to center the Abbe condenser in the optical light path.

- a) Close the field diaphragm (7). If the condenser is out of center, you will see an image of the iris opening similar to Fig. A. The image may be out of focus and out of center.
- b) Adjust the condenser up-down knob (1) until the side of image is clear. In order to center the iris opening in the field of view, use the two centering screws (5) at each side of the Abbe condenser mounting ring (see Fig. B).
- c) Once the image is in focus and centered in the field of view, open the field diaphragm (7) until the image is at least as large as the field of view (see Fig. C).

- 3) The field diaphragm (7) on the lamp house has a filter seat (8) for holding the filters. You may put a filter on it to enhance the performance.
- 4) When you watch with the 4x objective, put the auxiliary lens of condenser (2) into observation position. You may adjust the aperture diaphragm (with 3) when you use different objectives.

#### 6. Polarization observation (Fig. 8)

1. The polarizing observation function is selected by pushing the analyzer (5) into the light path and using the polarizer (Fig. 7, No. 6). Rotate both of them to „0“ (e. g. the analyzer with 2); that may let the polarizer and analyzer be orthogonally, namely the polarizer with E-W direction of passage, the analyzer with N-S direction of passage.

8

- 1 Rod for observation/photo changing
- 2 Analyzer adjustment tray
- 3 Bertrand lens centering screw
- 4 Compensator
- 5 Analyzer


2. To make path differences, you may push the compensator,  $\lambda/4$  compensator or quartz wedge compensator (4) into the light path.
3. The push rod to the upper left of the trinocular head is the rod for observation/photo changing. In the observation mode, 100% of the image light will be sent to the binocular viewing; in the photo mode, the light will be sent to the camera port on the top. It satisfies high-grape photography.

**Note:** If the condenser is adjusted, you will adjust the polarizer again (see Fig. 8). Set the analyzer to light path and rotate to „0“. Rotate the polarizer to „0“. Loosen the condenser holding screw, rotate the polarizer mark until the field is very dark. Then lock the screw again.

### 7. Conoscopic observation (Fig. 8)

Conoscopic observation is selected by pushing the Bertrand lens (Fig. 2, No. 3) into the light path. Adjusting the Bertrand lens centering screw (3) may let the Bertrand lens in the center of view.

### 8. Exchange of the lamp (Fig. 9)


- 1 Lamp socket
- 2 Bulb
- 3 Lamp door retaining screw
- 4 Lamp centering screw

1. Disconnect the power plug and wait to be sure the bulb is cool. Open the lamp door using the lamp door retaining screw (3).

2. Pull the old bulb (2) straight out of the socket (1). Insert the new one into the same fixture. When you install the new lamp, be careful not to touch the glass with your fingers. The new bulb should be supplied in a plastic protective envelope. If not, use a tissue or another medium to grasp the bulb. This will prevent contamination from your hand, which can reduce the lamp's intensity and life. Reinsert the lamp socket into the lamp house and retighten the screw (3). If necessary, you may need to adjust the lamp centering screw (4) to center the bulb.

### 9. Exchange of the fuse (Fig. 10)

Switch off the power switch and pull out the plugs of the power cord from mains socket and from mains in at the microscope (1). Grip with a small screwdriver under the lug of fuse holder (2) and pull it carefully forwards (Fig. 10, arrow). Take out the fuse holder

- (2). Change the fuse. Install the fuse holder and plug again. The specification of fuse: 250 V / 0,5 A, F.

## V. MAINTENANCE

### 1. Cleaning frame and stage

Disconnect the plug from mains socket before cleaning. Clean the frame and the stage with a soft cloth moistened with a mild detergent solution. Be sure that the instrument is dry before you use it.

### 2. Cleaning optical parts

Microscope eyepieces and objectives are coated. They should not be wiped while dry dirt or dust may scratch the coating. It is best to remove parts from the frame prior to cleaning. Always blow loose dust away first. Use lens tissue of good quality moistened with a lens cleaner or a small amount of alcohol, then wipe the surface clean with a lens tissue. Solvents such as xylene should NOT be used as cleaner.

### 3. Cleaning 100x oil immersion lens

The immersion oil should be removed from the lens at the end of each workday using lens tissue moistened with a lens cleaner or a small amount of alcohol.

**Do not disassemble objective lenses!**

**Bresser GmbH**

Gutenbergstr. 2 · DE-46414 Rhede · Germany

Tel. +49 (0) 28 72 - 80 74-210

Fax +49 (0) 28 72 - 80 74-222

[www.bresser.de](http://www.bresser.de) · [service@bresser.de](mailto:service@bresser.de)