



RED39-Z

No other stereo microscope stand can be compared either in design or quality. The base is completely sealed to protect the illuminator from liquids that might be spilled. All the models are fully equipped with incident and transmitted LED illumination, carrying handle for easy transportation and "one touch" sample clips.

Model RED 39-Z

Optical System Greenough

Observation Tube Binocular head

Interpupillary distance 51-77mm

Diopter adjustment On both tubes, +/- 5 diopter

Inclination 45° inclined, 360° rotating

Eyepieces WF10X/20mm

Objectives system Zoom, ratio 3:1

Objectives magnification 1X - 3X

Working distance 60mm

Stand Fixed arm type stand

Base 225x179mm

Column/Arm 202mm high

Head holder For Ø66mm stereo head with LED incident illumination

Focus mechanism Coarse focusing system with tension adjustment

Focusing stoke 46mm

Incident illuminationLED 3.5V/100mW with intensity controlTransmitted illuminationLED 3.5V/100mW with intensity control

Transformer External

Power supply Charger 100-240V or rechargeable batteries

Accessories included Black, white and frosted stage plates, dust cover, charger and batteries

Length x Width x Height 225x179x286mm

Weight 3.5Kg



Add a Moticam camera to capture, document, annotate and share images and videos with the Motic Images Plus 2.0 software.